

SUBJECT INDEX*

CONFORMATIONS

- 318, D320, 344, 347, 349, D357, 362, 363
- temp.-dependent, 283, 292, D293, D294

DISSOCIATION CONSTANTS, K_s

- substrates with
 - — microsomal P450, T225, 225
 - — solubilized P450, T265

ELECTRON

- carrier, 183, 196, D205, 223, 305, 373, 377
- transfer, F310, F311, 311, 313, 315, 319

ENTHALPY-ENTROPY

- (of) activation, T226, 290
- — compensation, 227, 228, 229, 282
- (of) substrate binding, 267, T268

HEME

- autooxidation, 295ff.
- — kinetics, F298, T299
- — spectra, solvent effects, F300, F301, 304
- axial ligands
 - — alkyl isocyanides, 282, F283, 340, F342, T343
 - — cyanide, 252, F253, 256, D260, T268, 340
 - — cysteYL, 244, 253
 - — energy diagrams, F241, F243
 - — ESR evidence, 252, F253
 - — H₂O, D249, 328, D329, D330, 219 (—OH)
 - — imidazole, 256, 267, T268, 340
 - — mercaptan (sulfur), F243, 244, 246, F247, D249, D250
 - — mercaptide (sulfur) 219, F243, 244, 246, F247, D249, D250, 337
 - — native, 194, 256, 267, 283, D320, 328, 337, 344
 - — pyridine, 198
 - — replacement, 189, T190, 259
- electronic structure, D249
- iron accessibility, 272, 290, D293, D294, 328, 336, 344, 345
- orientation, 239
- reduction kinetics, 301, T302

HEMOPROTEINS (COMPARISON)

- T185, F242, 245, D248, D250, 291, D293, T336, 337, T345

HYDROPHOBIC INTERACTIONS

- 267, D270, 273, 284, 318, 344, 353, 357, D357, D387

(HYDRO)PEROXIDES

- heme oxidation, 296, 302, 303
- P450 substrate oxidation, 168, 169, D175, D176, 200
- — kinetics, T228, D235

INDUCING AGENTS

- D269, 284, D364

INTERACTIONS

- solute-solute, 311, 315, 318, D320, 347ff., 353ff., F355, D357
- solute-solvent, 347, 354

KINETICS, TEMP. DEPENDENCE

- T226

MECHANISM (HYDROXYLATION)

- 172, 183, 221
- intermediates, 170, F170
- oxenoid, 220, 230, 231, 232, 233
- radicals, 233, D235
- sequential, F275
- spin equilibrium, 259

MEMBRANE(S)

- F349, 360, 363, D365, D366, 370, 373
- enzyme activities 374, T376, 377
- enzyme complexes, 385
- intramembrane interactions
 - — accross-membrane, F386
 - — attractive/repulsive, 353, 381, F385
 - — elastic, 350, 380, F381, 386, D387, D388
 - — long/short range, 353, 384, D387
 - — non-Van der Waals, 382
 - — protein-lipid, F349, 373, 379ff., 384, D387
 - — protein-protein, T350, 379ff., 384
- phase transitions, 348, T350, D357, 385, D387
- reconstitution, 368ff., F369, F370, F371, T372, T374
- solubilization, F368

METHODS

- absorption spectroscopy, 181, 264, F265

* D = discussion, F = figure, T = table

X

- — polarized, 238
- — temp.-differ., 252, 256
- circular dichroism 339ff.
- ESR, 181, 187, 252
- fast kinetics, 272
- Hückel, extended, 238
- infrared spectroscopy, 334
- preparative (biochemical), 164, 166, 181, 184, 185, 207ff., 359, 360, 368
- proton magnetic relaxation, 324
- stopped-flow kinetics, 297
- sub-zero temperature (solvents)
- — (in) absorption spectroscopy, 279, 310
- — (in) photodissoc. kinetics, 287

MICROSOMES

- 360, F361, T361, 363, 365

MODELS

- hydroxylation
- — metal complexes +O₂, 231, T231
- — peracids of metals (no O₂), 231, T232
- structural (active site), F216, F217, F218, T218, D294, 296, 297, 299, 337, D338

MULTI-ENZYME SYSTEMS

- T169, 273, 319, 385

P420

- (in) P450ADR, 186
- (from) P450LM
- — absorp. spectra, F335
- — (by) phospholipase C treatment, F362, D366
- — (by) temp. treatment, F334, T336
- (by) protonation of cys in P450, D249
- reconversion into P450, F369, F370, F371, F373, 377

P450ADR

- 179ff.
- +CO (affinity), F197, D203
- +hydroperoxides, F200, F201
- +NO, 199
- +O₂, 195, 196
- perturbants, 191—193
- preparation, 184ff., D205
- properties, 185
- — ESR spectra, T185, T189
- — redox state, 194, D205
- — spin state, 186
- +substrate/ligand, F187, T189, F191, F198

P450CAM

- absorption spectra, polarized, F239, F240, D248
- axial ligands, T218, F242, 244, 245, 246, F247, D248, D249

- +CO F242, F243, D250
- dienzyme complex, 315, 319, D320
- +diethylphosphine, F215
- low-temp. studies, F280, F289, F290, 292
- +O₂, 220
- redox reactions, F310, F311, F314, D320
- spin equilibria, T281, F316, D320

P450LM

- +CO, F334 (abs. sp.); 340, F342 (CD); F334, T336 (IR); 272, F272 (kinet.); F325 (PMR)
- +ethylisocyanide, 283, F283, 340, F342, T343
- PMR-rates, F325
- reduction kinetics, T274
- subspecies, D260, 266, D269, 283, 329, 335, 361
- — characterization, 164, F165, 167, T168 (LM2); 164, F165, D175 (LM4)
- +substrates
- — binding kinetics, T275, D276
- — spectra, F254, T255, F256, F257, F258 (temp. diff.); F341, T343 (CD)

PHOSPHOLIPASE C

- 360ff., D357

PHOSPHO

LIPID(S)

- boundary (annular), 349, 363, D365, D387
- D203, D269, 273, 347, 348, F349, F352, D357, 360, F361, T361, 363, D365, D366, 368, 371, T372, 373, 379ff., F384, 385, D387

PROTONATION

- 220, 244, 245
- spin equil., 282, D285

PUTIDAREDOXIN(Pd)-P450CAM(m) COMPLEX

- 310, 311; 313 (assoc. const. Pd^o-m^{os}; E_a); F314, 315 (decomp. rate); F317 (diss. const. Pd^o-m^{os}); 319, D320

RATE LIMITING STEP

- 229 (O₂ → substrate); 277 (substrate binding)

RELAXATION TIMES

- electronic, 252, 328, D330
- PMR, D294, 323ff., 327 (correl. time)
- spin equil., 259

SOLVENT

- correlation length, 347, 350, 352, 353
- order parameter, 348, F349, F352, 354

- phase separation, 350
- structure (perturbation) 347ff., 353ff.

SPIN STATE

- 165, 197, 244 (Fe^{3+} , 5/2); 165, 194, 219, 266, D320, 345 (Fe^{3+} , 1/2); 194, 244, D346 (Fe^{2+} , 4/2); 245, 246 (Fe^{2+} , 0)
- equilibria, 252, 254, F256, F257, F258, 259, 280ff., D285, 315, D320
- kinetics, F316, D320
- pH dependence, F193
- thermodynamics, 258, T281

SUBSTRATE BINDING

- 344
- aniline, F273
- deoxycorticosterone, 187, F188
- ESR spectra, 252, F253
- kinetics, 271
- progesterone, F209
- temp. diff. spectra, F256, F257, F258
- thermodynamics, 266, T267
- type I/type II, 258, 266, F267, 267, T268, D269, D270, T274, 340, F341, T343